





Cryometrix T-160 Ultra Low Temp Freezer

Download Brochure

Calculate Storage Capacity

Download Map

The new 2017 redesigned Cryometrix T-160 programmable ULT freezer delivers on all points. The T-160 freezer is a technological breakthrough that incorporates a compressor-free, redundant, liquid nitrogen-based technology. The T-160 is a stable, scalable, adaptable, and validatable platform. The payload capacity is 33.5 ft3 (950 L = 72,000 2-mL vials), and with its all stainless steel upright design, the T-160 addresses several user issues such as security, ease of operation, cleaning, and the system is designed to eliminate icing.

Physical Property	Units	Values
Dimensions	Inches	42" D x 48" W x 84" H
Weight (Empty)	Pounds	1672
Payload Capacity	Cubic Feet	31.7
Payload Dimensions	Inches	26.5" D x 34.5" W x 60" H
Electrical	120 VAC 60 Hz	
Cooling Media	LN2	18-90
Inlet Pressure	PSI	+68 to -256
Adjustable Set Point	°F	± 7°F at -160°C
Temperature Uniformity	°F	± 5°F at -80°C
Ramp Rate	°F/min	-2 to -5

- Fast controllable temperature ramp

Unique sealed Liquid Nitrogen based technology

- Uniform temperatures throughout payload bay (± 5°C)
- Fast recovery from open/close event (10 mins)
- Large 31.7 cubic ft. / 900 liters storage volume

- 57,600 x 2ml vial storage capacity controlled at \pm 3°C

- Customized racking available
- Compressor free operation

- Minimal power consumption - No rejected heat into workspace

- Noiseless operation, vibration and dust free

Low LN2 consumption Mission critical design - Redundant circuits assure reliability

- User security features - Low maintenance design
- Ergonomically simplified design - As easy to use as a common household refrigerator

- Left or right hand door swing

- Designed for cGMP applications - DQ/IQ completed
- OQ/PQ ready - GMP qualified units currently in PharmBio use
- Compared to compressor based technologies The T-160 offers:

- A significant reduction in electrical power consumption

- Up to 5 times greater payload capacity per unit area of floor space - Significant overall reduction in net operating expenses
- No liquid nitrogen exposure to users or product

Products T-160 Ultra Low Freezer

T-90 Storage Freezer

S-90 Shipper Freezer

B-90 Blast/Thaw Freezer

Not seeing what you're looking for? Need something custom? Please reach out, and we'll see what we can do!

Send us a Message



Cryometrix T-160 Ultra Low Temp Freezer

<u>Download Brochure</u> <u>Calculate Storage Capacity</u> <u>Download Map</u>

The new 2017 redesigned Cryometrix T-160 programmable ULT freezer delivers on all points. The T-160 freezer is a technological breakthrough that incorporates a compressor-free, redundant, liquid nitrogen-based technology. The T-160 is a stable, scalable, adaptable, and validatable platform. The payload capacity is 33.5 ft3 (950 L = 72,000 2-mL vials), and with its all stainless steel upright design, the T-160 addresses several user issues such as security, ease of operation, cleaning, and the system is designed to eliminate icing.

Physical Property	Units	Values
Dimensions	Inches	42" D x 48" W x 84" H
Weight (Empty)	Pounds	1672
Payload Capacity	Cubic Feet	31.7
Payload Dimensions	Inches	26.5" D x 34.5" W x 60" H
Electrical	120 VAC 60 Hz	
Cooling Media	LN2	18-90
Inlet Pressure	PSI	+68 to -256
Adjustable Set Point	°F	± 7°F at -160°C
Temperature Uniformity	°F	± 5°F at -80°C
Ramp Rate	°F/min	-2 to -5



Features & Benefits

Unique sealed Liquid Nitrogen based technology

- Fast controllable temperature ramp
- Uniform temperatures throughout payload bay (± 5°C)
- Fast recovery from open/close event (10 mins)

Large 31.7 cubic ft. / 900 liters storage volume

- 57,600 x 2ml vial storage capacity controlled at ± 3°C
- Customized racking available

Compressor free operation

- Minimal power consumption
- No rejected heat into workspace
- Noiseless operation, vibration and dust free

Low LN2 consumption Mission critical design

- Redundant circuits assure reliability
- User security features
- Low maintenance design

Ergonomically simplified design

- As easy to use as a common household refrigerator
- Left or right hand door swing

Designed for cGMP applications

- DQ/IQ completed
- OQ/PQ ready
- GMP qualified units currently in PharmBio use

Compared to compressor based technologies - The T-160 offers:

- A significant reduction in electrical power consumption
- Up to 5 times greater payload capacity per unit area of floor space
- Significant overall reduction in net operating expenses

No liquid nitrogen exposure to users or product

Products

T-160 Ultra Low Freezer
B-90 Blast/Thaw Freezer
T-90 Storage Freezer
S-90 Shipper Freezer

Not seeing what you're looking for? Need something custom? Please reach out, and we'll see what we can do!

Send us a Message